



EEG graph

Before High Data




Before Low Data




After High Data

Total Segments: 158
Hi Count: 101 (63.92%)
Lo Count: 57 (36.08%)



After Low Data

Total Segments: 158
Hi Count: 0 (0.00%)
Lo Count: 158 (100.00%)





RENDERED GRAPHS:

THE SYSTEM VISUALIZES THE RAW DATA OF THE SELECTED PARTICIPANT, DISPLAYING HIGH (HI) AND LOW (LOW) VALUES IN A GRAPH. ADDITIONALLY, IT RENDERS THE PREPROCESSED DATA AND MODEL PREDICTIONS FOR HIGH (HI) AND LOW (LOW) VALUES.

INTERACTIVE GRAPH FEATURES:

THE GRAPHS CAN BE ZOOMED IN AND OUT USING A MOUSE FOR BETTER ANALYSIS. TOOLTIPS ARE PROVIDED FOR EACH GRAPH SEGMENT, OFFERING DETAILED INFORMATION ABOUT THE DATA POINTS.

SHAP Values (Channel-wise)




VISUALIZATION OF GRAPHS AND NORMALIZED SHAP VALUES

1. GRAPH VISUALIZATION:

- THE NORMALIZED SHAP VALUES FOR EACH PREPROCESSED CHANNEL ARE VISUALIZED.

2. LABEL INFORMATION:

- CNN VERSION: DISPLAYS VALUES FOR A TOTAL OF 14 LABELS.
- CNN_FFT VERSION: DISPLAYS VALUES FOR A TOTAL OF 70 LABELS, CORRESPONDING TO CHANNEL + BANDWIDTH COMBINATIONS.



Contour Visualization Based on Error Range

Each channel's contour map is displayed using a red color scheme, where the larger the error range between the maximum (max) and minimum (min) values of the EEG data, the more intense the red color appears.